

## Synaptotagmin4

Cat.No. 105 143; Polyclonal rabbit antibody, 50 µg specific antibody (lyophilized)

### Data Sheet

Reconstitution/ Storage	50 µg specific antibody, lyophilized. Affinity purified with the immunogen. Albumin and azide were added for stabilization. For <b>reconstitution</b> add 50 µl H <sub>2</sub> O to get a 1mg/ml solution in PBS. Then aliquot and store at -20°C to -80°C until use. Antibodies should be stored at +4°C when still lyophilized. Do not freeze! For detailed information, see back of the data sheet.
Applications	<b>WB:</b> 1 : 1000 (AP staining) (see remarks) <b>IP:</b> not tested yet <b>ICC:</b> 1 : 500 up to 1 : 2000 <b>IHC:</b> 1 : 200 up to 1 : 500 <b>IHC-P:</b> not tested yet
Immunogen	Recombinant protein corresponding to AA 40 to 151 from rat Synaptotagmin4 (UniProt Id: P50232)
Reactivity	Reacts with: human (Q9H2B2), rat (P50232), mouse (P40749), hamster. Other species not tested yet.
Specificity	K.O. validated
Matching control	105-4P
Remarks	<b>WB:</b> Cat. no. <a href="#">105 043</a> is recommended for this application.

**TO BE USED IN VITRO / FOR RESEARCH ONLY**  
**NOT TOXIC, NOT HAZARDOUS, NOT INFECTIOUS, NOT CONTAGIOUS**

### Background

Up to now at least 17 synaptotagmins have been identified. **Synaptotagmin4** is composed of a vesicular, a transmembrane and two C2 domains. Only the C2B domain is able to bind calcium. In the C2A domain one of the calcium binding aspartates has been substituted for serine leading to a loss of its binding capabilities.

The localization of synaptotagmin4 is still under discussion. A localization to synaptic vesicles (SVs) has been postulated but more recent studies suggest that it is present in the Golgi compartment, in distal parts of neurites and on large dense core vesicles (LDCVs) of NGF differentiated PC12 cells.

### Selected References for 105 143

A novel method for culturing stellate astrocytes reveals spatially distinct Ca<sup>2+</sup> signaling and vesicle recycling in astrocytic processes.

Wolfe AC, Ahmed S, Awasthi A, Stahlberg MA, Rajput A, Magruder DS, Bonn S, Dean C  
The Journal of general physiology (2017) 149: 149-170. . **ICC, WB**

Conserved expression of the zebrafish syt4 gene in GABAergic neurons in the cerebellum of adult fishes revealed by mammalian SYT4 immunoreactive-like signals.

Shiao MS, Liu ST, Siriwatcharapibool G, Thongpradit S, Khunpanich P, Tong SK, Huang CH, Jinawath N, Chou MY  
Heliyon (2024) 10: e30575. . **IHC; tested species: zebrafish**

Sex-specific regulation of follicle-stimulating hormone secretion by synaptotagmin 9.

Roper LK, Briguglio JS, Evans CS, Jackson MB, Chapman ER  
Nature communications (2015) 6: 8645. . **ICC**

### Selected General References

Altered hippocampal short-term plasticity and associative memory in synaptotagmin IV (-/-) mice.  
Ferguson GD et al. Hippocampus (2004) PubMed:15390175

Structural basis for the evolutionary inactivation of Ca<sup>2+</sup> binding to synaptotagmin 4.  
Dai H et al. Nat. Struct. Mol. Biol. (2004) PubMed:15311271

Reduced anxiety and depression-like behavior in synaptotagmin IV (-/-) mice.  
Ferguson GD et al. Neuropharmacology (2004) PubMed:15380377

Synaptotagmin IV regulates glial glutamate release.  
Zhang Q et al. Proc. Natl. Acad. Sci. U.S.A. (2004) PubMed:15197251

Nerve growth factor-dependent sorting of synaptotagmin IV protein to mature dense-core vesicles that undergo calcium-dependent exocytosis in PC12 cells.  
Fukuda M et al. J. Biol. Chem. (2003) PubMed:12446703

Non-polarized distribution of synaptotagmin IV in neurons: evidence that synaptotagmin IV is not a synaptic vesicle protein.  
Ibata K et al. Neurosci. Res. (2002) PubMed:12135783

A unique spacer domain of synaptotagmin IV is essential for Golgi localization.  
Fukuda M et al. J. Neurochem. (2001) PubMed:11331402

Synaptotagmin IV is present at the Golgi and distal parts of neurites.  
Ibata K et al. J. Neurochem. (2000) PubMed:10646502

Functional and biochemical analysis of the C2 domains of synaptotagmin IV.  
Thomas DM et al. Mol. Biol. Cell (1999) PubMed:10397765

Synaptotagmin IV: biochemistry, genetics, behavior, and possible links to human psychiatric disease.  
Ferguson GD et al. Mol. Neurobiol. () PubMed:11817218

Access the online factsheet including applicable protocols at <https://sysy.com/product/105143> or scan the QR-code.



## FAQ - How should I store my antibody?

### Shipping Conditions

- All our antibodies and control proteins / peptides are shipped lyophilized (vacuum freeze-dried) and are stable in this form without loss of quality at ambient temperatures for several weeks.

### Storage of Sealed Vials after Delivery

- **Unlabeled** and **biotin-labeled antibodies** and **control proteins** should be stored at 4°C before reconstitution. **They must not be stored in the freezer when still lyophilized!** Temperatures below zero may cause loss of performance.
- **Fluorescence-labeled antibodies** should be reconstituted immediately upon receipt. Long term storage (several months) may lead to aggregation.
- **Control peptides** should be kept at -20°C before reconstitution.

### Long Term Storage after Reconstitution (General Considerations)

- The storage freezer must not be of the frost-free variety ("no-frost freezer"). This cycle between freezing and thawing (to reduce frost-build-up), which is exactly what should be avoided. For the same reason, antibody vials should be placed in an area of the freezer that has minimal temperature fluctuations, for instance towards the back rather than on a door shelf.
- Aliquot the antibody and store frozen (-20°C to -80°C). Avoid very small aliquots (below 20 µl) and use the smallest storage vial or tube possible. The smaller the aliquot, the more the stock concentration is affected by evaporation and adsorption of the antibody to the surface of the storage vial or tube. Adsorption of the antibody to the surface leads to a substantial loss of activity.
- The addition of glycerol to a final concentration of 50% lowers the freezing point of your stock and keeps your antibody at -20°C in liquid state. This efficiently avoids freeze and thaw cycles.

### Product Specific Hints for Storage

#### Control proteins / peptides

- Store at -20°C to -80°C.

#### Monoclonal Antibodies

- **Ascites** and **hybridoma supernatant** should be stored at -20°C up to -80°C. **Prolonged storage at 4°C is not recommended!** Unlike serum, ascites may contain proteases that will degrade the antibodies.
- **Purified IgG** should be stored at -20°C up to -80°C. Adding a carrier protein like BSA will increase long term stability. Many of our antibodies already contain carrier proteins. Please refer to the data-sheet for detailed information.

#### Polyclonal Antibodies

- **Crude antisera:** With anti-microbials added, they may be stored at 4°C. However, frozen storage (-20°C up to -80°C) is preferable.
- **Affinity purified antibodies:** Less robust than antisera. Storage at -20°C up to -80°C is recommended. Adding a carrier protein like BSA will increase long term stability. Most of our antibodies already contain carrier proteins. Please refer to the data-sheet for detailed information.

#### Fluorescence-labeled Antibodies

- Store as a liquid with 1 : 1 (v/v) glycerol at -20°C. Protect these antibodies from light exposure.

## Avoid repeated freeze-thaw cycles for all antibodies!

## FAQ - How should I reconstitute my antibody?

### Reconstitution

- All our purified antibodies are lyophilized from PBS. To reconstitute the antibody in PBS, add the amount of deionized water given in the respective datasheet. If higher volumes are preferred, add water as mentioned above and then the desired amount of PBS and a stabilizing carrier protein (e.g. BSA) to a final concentration of 2%. Some of our antibodies already contain albumin. Take this into account when adding more carrier protein. For complete reconstitution, carefully remove the lid. After adding water, briefly vortex the solution. You can spin down the liquid by placing the vial into a 50 ml centrifugation tube filled with paper.
- If desired, add small amounts of azide or thimerosal to prevent microbial growth. This is especially recommended if you want to keep an aliquot at 4°C.
- After reconstitution of fluorescence-labeled antibodies, add 1 : 1 (v/v) glycerol to a final concentration of 50%. This lowers the freezing point of your stock and keeps your antibody in liquid state at -20°C.
- Glycerol may also be added to unlabeled primary antibodies. It is a suitable way to avoid freeze-thaw cycles.
- Please refer to our **tips and hints for subsequent storage** of reconstituted antibodies and control peptides and proteins.