

## Mint1

Cat.No. 144 103; Polyclonal rabbit antibody, 50 µg specific antibody (lyophilized)

### Data Sheet

Reconstitution/ Storage	50 µg specific antibody, lyophilized. Affinity purified with the immunogen. Albumin and azide were added for stabilization. For <b>reconstitution</b> add 50 µl H <sub>2</sub> O to get a 1mg/ml solution in PBS. Then aliquot and store at -20°C to -80°C until use. Antibodies should be stored at +4°C when still lyophilized. Do not freeze! For detailed information, see back of the data sheet.
Applications	<b>WB:</b> 1 : 1000 (AP staining) <b>IP:</b> not tested yet <b>ICC:</b> 1 : 200 up to 1 : 500 <b>IHC:</b> 1 : 200 <b>IHC-P:</b> not tested yet
Immunogen	Recombinant protein corresponding to AA 2 to 265 from rat Mint1 (UniProt Id: O35430)
Reactivity	Reacts with: rat (O35430), mouse (B2RUJ5). Other species not tested yet.
Matching control	144-1P

**TO BE USED IN VITRO / FOR RESEARCH ONLY**  
**NOT TOXIC, NOT HAZARDOUS, NOT INFECTIOUS, NOT CONTAGIOUS**

## Background

Mints (also referred to as X11-like proteins) are adaptor-proteins that consist of divergent N-terminal sequences and conserved C-terminal PTB and PDZ domains.

Three isoforms (**Mint 1, 2 and 3**) have been described. Mint 1 exclusively binds to the synaptic protein CASK via its N-terminal sequence. Munc18-1 has been shown to be an interaction partner of Mint 1 and 2. The more C-terminal located PTB and PDZ domains present in all Mint isoforms bind to widely distributed proteins like APP, presenilins and Ca<sup>2+</sup> channels.

A Mint 1 knock out had no obvious effect on brain architecture and development, nor was synaptic plasticity in excitatory synapses affected. In inhibitory synapses of knock out strains the release of gamma-aminobutyric acid (GABA) was impaired.

## Selected References for 144 103

Molecular profiling of synaptic vesicle docking sites reveals novel proteins but few differences between glutamatergic and GABAergic synapses.

Boyken J, Grønborg M, Riedel D, Urlaub H, Jahn R, Chua JJ  
Neuron (2013) 782: 285-97. . **WB**

Liprin-α proteins are master regulators of human presynapse assembly.  
Marcó de la Cruz B, Campos J, Molinaro A, Xie X, Jin G, Wei Z, Acuna C, Sterky FH  
Nature neuroscience (2024) : . . **WB; tested species: human**

BACE1 controls synaptic function through modulating release of synaptic vesicles.  
Das B, Singh N, Yao AY, Zhou J, He W, Hu X, Yan R  
Molecular psychiatry (2021) 2611: 6394-6410. . **WB; tested species: mouse**

A role for Sec8 in oligodendrocyte morphological differentiation.  
Anitei M, Ifrim M, Ewart MA, Cowan AE, Carson JH, Bansal R, Pfeiffer SE  
Journal of cell science (2006) 119Pt 5: 807-18. . **WB; tested species: rat**

## Selected General References

Mint1, a Munc-18-interacting protein, is expressed in insulin-secreting beta-cells.  
Zhang W, Lilja L, Bark C, Berggren PO, Meister B  
Biochemical and biophysical research communications (2004) 3203: 717-21. .

A multiprotein trafficking complex composed of SAP97, CASK, Veli, and Mint1 is associated with inward rectifier Kir2 potassium channels.

Leonoudakis D, Conti LR, Radeke CM, McGuire LM, Vandenberg CA  
The Journal of biological chemistry (2004) 27918: 19051-63. .

A role for Mints in transmitter release: Mint 1 knockout mice exhibit impaired GABAergic synaptic transmission.  
Ho A, Morishita W, Hammer RE, Malenka RC, Sudhof TC  
Proceedings of the National Academy of Sciences of the United States of America (2003) 1003: 1409-14. .

Regulation of APP-dependent transcription complexes by Mint/X11: differential functions of Mint isoforms.  
Biederer T, Cao X, Südhof TC, Liu X  
The Journal of neuroscience : the official journal of the Society for Neuroscience (2002) 2217: 7340-51. .

CASK participates in alternative tripartite complexes in which Mint 1 competes for binding with caskin 1, a novel CASK-binding protein.

Tabuchi K, Biederer T, Butz S, Sudhof TC  
The Journal of neuroscience : the official journal of the Society for Neuroscience (2002) 2211: 4264-73. .

Amyloid precursor protein associates independently and collaboratively with PTB and PDZ domains of mint on vesicles and at cell membrane.

Okamoto M, Nakajima Y, Matsuyama T, Sugita M  
Neuroscience (2001) 1043: 653-65. .

Access the online factsheet including applicable protocols at <https://sysy.com/product/144103> or scan the QR-code.



# FAQ - How should I store my antibody?

## Shipping Conditions

- All our antibodies and control proteins / peptides are shipped lyophilized (vacuum freeze-dried) and are stable in this form without loss of quality at ambient temperatures for several weeks.

## Storage of Sealed Vials after Delivery

- **Unlabeled** and **biotin-labeled antibodies** and **control proteins** should be stored at 4°C before reconstitution. **They must not be stored in the freezer when still lyophilized!** Temperatures below zero may cause loss of performance.
- **Fluorescence-labeled antibodies** should be reconstituted immediately upon receipt. Long term storage (several months) may lead to aggregation.
- **Control peptides** should be kept at -20°C before reconstitution.

## Long Term Storage after Reconstitution (General Considerations)

- The storage freezer must not be of the frost-free variety ("no-frost freezer"). This cycle between freezing and thawing (to reduce frost-build-up), which is exactly what should be avoided. For the same reason, antibody vials should be placed in an area of the freezer that has minimal temperature fluctuations, for instance towards the back rather than on a door shelf.
- Aliquot the antibody and store frozen (-20°C to -80°C). Avoid very small aliquots (below 20 µl) and use the smallest storage vial or tube possible. The smaller the aliquot, the more the stock concentration is affected by evaporation and adsorption of the antibody to the surface of the storage vial or tube. Adsorption of the antibody to the surface leads to a substantial loss of activity.
- The addition of glycerol to a final concentration of 50% lowers the freezing point of your stock and keeps your antibody at -20°C in liquid state. This efficiently avoids freeze and thaw cycles.

## Product Specific Hints for Storage

### Control proteins / peptides

- Store at -20°C to -80°C.

### Monoclonal Antibodies

- **Ascites** and **hybridoma supernatant** should be stored at -20°C up to -80°C. **Prolonged storage at 4°C is not recommended!** Unlike serum, ascites may contain proteases that will degrade the antibodies.
- **Purified IgG** should be stored at -20°C up to -80°C. Adding a carrier protein like BSA will increase long term stability. Many of our antibodies already contain carrier proteins. Please refer to the data-sheet for detailed information.

### Polyclonal Antibodies

- **Crude antisera:** With anti-microbials added, they may be stored at 4°C. However, frozen storage (-20°C up to -80°C) is preferable.
- **Affinity purified antibodies:** Less robust than antisera. Storage at -20°C up to -80°C is recommended. Adding a carrier protein like BSA will increase long term stability. Most of our antibodies already contain carrier proteins. Please refer to the data-sheet for detailed information.

### Fluorescence-labeled Antibodies

- Store as a liquid with 1 : 1 (v/v) glycerol at -20°C. Protect these antibodies from light exposure.

# Avoid repeated freeze-thaw cycles for all antibodies!

## FAQ - How should I reconstitute my antibody?

### Reconstitution

- All our purified antibodies are lyophilized from PBS. To reconstitute the antibody in PBS, add the amount of deionized water given in the respective datasheet. If higher volumes are preferred, add water as mentioned above and then the desired amount of PBS and a stabilizing carrier protein (e.g. BSA) to a final concentration of 2%. Some of our antibodies already contain albumin. Take this into account when adding more carrier protein. For complete reconstitution, carefully remove the lid. After adding water, briefly vortex the solution. You can spin down the liquid by placing the vial into a 50 ml centrifugation tube filled with paper.
- If desired, add small amounts of azide or thimerosal to prevent microbial growth. This is especially recommended if you want to keep an aliquot a 4°C.
- After reconstitution of fluorescence-labeled antibodies, add 1 : 1 (v/v) glycerol to a final concentration of 50%. This lowers the freezing point of your stock and keeps your antibody in liquid state at -20°C.
- Glycerol may also be added to unlabeled primary antibodies. It is a suitable way to avoid freeze-thaw cycles.
- Please refer to our **tips and hints for subsequent storage** of reconstituted antibodies and control peptides and proteins.