

GABA-A receptor $\alpha 2$ extracellular

Cat.No. 224 103; Polyclonal rabbit antibody, 50 μ g specific antibody (lyophilized)

Data Sheet

Reconstitution/ Storage	50 μ g specific antibody, lyophilized. Affinity purified with the immunogen. Albumin and azide were added for stabilization. For reconstitution add 50 μ l H ₂ O to get a 1mg/ml solution in PBS. Then aliquot and store at -20°C to -80°C until use. Antibodies should be stored at +4°C when still lyophilized. Do not freeze! For detailed information, see back of the data sheet.
Applications	WB: yes (see remarks) IP: yes ICC: 1 : 500 (see remarks) IHC: 1 : 1000 (see remarks) IHC-P: not tested yet EM: yes
Immunogen	Synthetic peptide corresponding to AA 29 to 37 from rat GABA-A receptor $\alpha 2$ (UniProt Id: P23576)
Reactivity	Reacts with: human (P47869), rat (P23576), mouse (P26048). Other species not tested yet.
Specificity	K.O. validated PubMed: 26973458
Matching control	224-1P
Remarks	WB: The crude antiserum (cat. no. 224 102) is highly recommended. GABA-A receptor alpha2 aggregates after boiling, making it necessary to run SDS-PAGE with non-boiled samples. ICC: This antibody is also suitable for the surface staining of living cells. After washing cells with bound antibodies, they can be fixed and visualized with secondary reagents. IHC: For best results use the protocol of Schneider Gasser et al., 2006.

TO BE USED IN VITRO / FOR RESEARCH ONLY
NOT TOXIC, NOT HAZARDOUS, NOT INFECTIOUS, NOT CONTAGIOUS

Background

Gamma-aminobutyric acid type **A (GABA-A)** receptors mediate the majority of inhibitory neurotransmission in the brain. These receptor proteins are ligand gated chloride ion channels and consist of a pentameric combination of different subunits (**alpha**, beta, gamma, delta, epsilon and rho). The resulting heterogenous population of GABA-A receptor subtypes are expressed throughout the brain with specific cellular and subcellular expression patterns.

Selected References for 224 103

- GABA(A) receptors containing the $\alpha 2$ subunit are critical for direction-selective inhibition in the retina. Auferkorte ON, Baden T, Kaushalya SK, Zabouri N, Rudolph U, Haverkamp S, Euler T PLoS one (2012) 74: e35109. . **IHC, EM; tested species: mouse, rabbit**
- DNA repair enzyme NEIL3 enables a stable neural representation of space by shaping transcription in hippocampal neurons. Kunath N, Bugaj AM, Bigonah P, Fernandez-Berrocal MS, Bjørås M, Ye J iScience (2021) 2412: 103470. . **WB, IHC, FR; tested species: mouse**
- Identification of a Core Amino Acid Motif within the α Subunit of GABAARs that Promotes Inhibitory Synaptogenesis and Resilience to Seizures. Nathanson AJ, Zhang Y, Smalley JL, Ollerhead TA, Rodriguez Santos MA, Andrews PM, Wobst HJ, Moore YE, Brandon NJ, Hines RM, Davies PA, et al. Cell reports (2019) 283: 670-681.e8. . **WB, ICC; tested species: mouse**
- Early postnatal GABAA receptor modulation reverses deficits in neuronal maturation in a conditional neurodevelopmental mouse model of DISC1. Saito A, Taniguchi Y, Rannals MD, Merfeld EB, Ballinger MD, Koga M, Ohtani Y, Gurley DA, Sedlak TW, Cross A, Moss SJ, et al. Molecular psychiatry (2016) 2110: 1449-59. . **ICC, IHC**
- Inhibitory synapse dysfunction and epileptic susceptibility associated with KIF2A deletion in cortical interneurons. Ruiz-Reig N, García-Sánchez D, Schakman O, Gailly P, Tissir F Frontiers in molecular neuroscience (2022) 15: 1110986. . **WB, IHC; tested species: mouse**
- The $\alpha 3$ subunit of GABAA receptors promotes formation of inhibitory synapses in the absence of collybistin. Wagner S, Lee C, Rojas L, Specht CG, Rhee J, Brose N, Papadopoulos T The Journal of biological chemistry () 296: 100709. . **ICC, IHC; tested species: mouse**
- Neural extracellular matrix regulates visual sensory motor integration. Reinhard J, Mueller-Buehl C, Wiemann S, Roll L, Luft V, Shabani H, Rathbun DL, Gan L, Kuo CC, Franzen J, Joachim SC, et al. iScience (2024) 272: 108846. . **IHC; tested species: mouse**
- Protocol for in vivo analysis of pre- and post-synaptic protein function in mice. Cramer TML, Abegg A, Tyagarajan SK STAR protocols (2024) 52: 103117. . **IHC; tested species: mouse**
- Sleep and wake cycles dynamically modulate hippocampal inhibitory synaptic plasticity. Wu K, Han W, Lu W PLoS biology (2022) 2011: e3001812. . **WB; tested species: mouse**
- REST/NRSF drives homeostatic plasticity of inhibitory synapses in a target-dependent fashion. Prestigio C, Ferrante D, Marte A, Romei A, Lignani G, Onofri F, Valente P, Benfenati F, Baldelli P eLife (2021) 10: . . **WB; tested species: mouse**
- NEIL1 and NEIL2 DNA glycosylases modulate anxiety and learning in a cooperative manner in mice. Hildrestrand GA, Rolseth V, Kunath N, Suganthan R, Jensen V, Bugaj AM, Fernandez-Berrocal MS, Sikko SB, Vetlesen S, Kuśnierczyk A, Olsen AK, et al. Communications biology (2021) 41: 1354. . **IHC; tested species: mouse**
- Tonic GABAergic inhibition, via GABAA receptors containing $\alpha \beta \epsilon$ subunits, regulates excitability of ventral tegmental area dopamine neurons. Tossell K, Dodhia RA, Galet B, Tkachuk O, Ungless MA The European journal of neuroscience (2021) . . **IHC; tested species: mouse**

Access the online factsheet including applicable protocols at <https://sysy.com/product/224103> or scan the QR-code.



FAQ - How should I store my antibody?

Shipping Conditions

- All our antibodies and control proteins / peptides are shipped lyophilized (vacuum freeze-dried) and are stable in this form without loss of quality at ambient temperatures for several weeks.

Storage of Sealed Vials after Delivery

- **Unlabeled** and **biotin-labeled antibodies** and **control proteins** should be stored at 4°C before reconstitution. **They must not be stored in the freezer when still lyophilized!** Temperatures below zero may cause loss of performance.
- **Fluorescence-labeled antibodies** should be reconstituted immediately upon receipt. Long term storage (several months) may lead to aggregation.
- **Control peptides** should be kept at -20°C before reconstitution.

Long Term Storage after Reconstitution (General Considerations)

- The storage freezer must not be of the frost-free variety ("no-frost freezer"). This cycle between freezing and thawing (to reduce frost-build-up), which is exactly what should be avoided. For the same reason, antibody vials should be placed in an area of the freezer that has minimal temperature fluctuations, for instance towards the back rather than on a door shelf.
- Aliquot the antibody and store frozen (-20°C to -80°C). Avoid very small aliquots (below 20 µl) and use the smallest storage vial or tube possible. The smaller the aliquot, the more the stock concentration is affected by evaporation and adsorption of the antibody to the surface of the storage vial or tube. Adsorption of the antibody to the surface leads to a substantial loss of activity.
- The addition of glycerol to a final concentration of 50% lowers the freezing point of your stock and keeps your antibody at -20°C in liquid state. This efficiently avoids freeze and thaw cycles.

Product Specific Hints for Storage

Control proteins / peptides

- Store at -20°C to -80°C.

Monoclonal Antibodies

- **Ascites** and **hybridoma supernatant** should be stored at -20°C up to -80°C. **Prolonged storage at 4°C is not recommended!** Unlike serum, ascites may contain proteases that will degrade the antibodies.
- **Purified IgG** should be stored at -20°C up to -80°C. Adding a carrier protein like BSA will increase long term stability. Many of our antibodies already contain carrier proteins. Please refer to the data-sheet for detailed information.

Polyclonal Antibodies

- **Crude antisera:** With anti-microbials added, they may be stored at 4°C. However, frozen storage (-20°C up to -80°C) is preferable.
- **Affinity purified antibodies:** Less robust than antisera. Storage at -20°C up to -80°C is recommended. Adding a carrier protein like BSA will increase long term stability. Most of our antibodies already contain carrier proteins. Please refer to the data-sheet for detailed information.

Fluorescence-labeled Antibodies

- Store as a liquid with 1 : 1 (v/v) glycerol at -20°C. Protect these antibodies from light exposure.

Avoid repeated freeze-thaw cycles for all antibodies!

FAQ - How should I reconstitute my antibody?

Reconstitution

- All our purified antibodies are lyophilized from PBS. To reconstitute the antibody in PBS, add the amount of deionized water given in the respective datasheet. If higher volumes are preferred, add water as mentioned above and then the desired amount of PBS and a stabilizing carrier protein (e.g. BSA) to a final concentration of 2%. Some of our antibodies already contain albumin. Take this into account when adding more carrier protein. For complete reconstitution, carefully remove the lid. After adding water, briefly vortex the solution. You can spin down the liquid by placing the vial into a 50 ml centrifugation tube filled with paper.
- If desired, add small amounts of azide or thimerosal to prevent microbial growth. This is especially recommended if you want to keep an aliquot a 4°C.
- After reconstitution of fluorescence-labeled antibodies, add 1 : 1 (v/v) glycerol to a final concentration of 50%. This lowers the freezing point of your stock and keeps your antibody in liquid state at -20°C.
- Glycerol may also be added to unlabeled primary antibodies. It is a suitable way to avoid freeze-thaw cycles.
- Please refer to our **tips and hints for subsequent storage** of reconstituted antibodies and control peptides and proteins.