

BASP1

Cat.No. 246 003; Polyclonal rabbit antibody, 50 µg specific antibody (lyophilized)

Data Sheet

Reconstitution/ Storage	50 µg specific antibody, lyophilized. Affinity purified with the immunogen. Albumin and azide were added for stabilization. For reconstitution add 50 µl H ₂ O to get a 1mg/ml solution in PBS. Then aliquot and store at -20°C to -80°C until use. Antibodies should be stored at +4°C when still lyophilized. Do not freeze! For detailed information, see back of the data sheet.
Applications	WB: 1 : 1000 (see remarks) IP: not tested yet ICC: 1 : 100 up to 1 : 500 IHC: not tested yet IHC-P: 1 : 500
Immunogen	Recombinant protein corresponding to AA 1 to 220 from rat BASP1 (UniProt Id: Q05175)
Reactivity	Reacts with: rat (Q05175), mouse. Other species not tested yet.
Remarks	WB: Due to its acidic pH, this protein migrates at higher molecular weight than theoretically calculated.

TO BE USED IN VITRO / FOR RESEARCH ONLY
NOT TOXIC, NOT HAZARDOUS, NOT INFECTIOUS, NOT CONTAGIOUS

Background

The **brain acid soluble protein 1 (BASP 1)** or **neuronal axonal membrane protein (NAP-22)** is a Ca²⁺-dependent calmodulin-binding protein that is specifically expressed in neuronal tissue. BASP 1 has several phosphorylation sites, undergoes N-terminal myristoylation and has been shown to be involved in intracellular signalling processes modulating neurite outgrowth.

Selected General References

Ion channel activity of brain abundant protein BASP1 in planar lipid bilayers.
Ostroumova OS, Schagina LV, Mosevitsky MI, Zakharov VV
The FEBS journal (2011) 2783: 461-9. .

Oligomeric structure of brain abundant proteins GAP-43 and BASP1.
Zakharov VV, Mosevitsky MI
Journal of structural biology (2010) 1703: 470-83. .

Dynamic interaction between WT1 and BASP1 in transcriptional regulation during differentiation.
Green LM, Wagner KJ, Campbell HA, Addison K, Roberts SG
Nucleic acids research (2009) 372: 431-40. .

Characterization of BASP1-mediated neurite outgrowth.
Korshunova I, Caroni P, Kolkova K, Berezin V, Bock E, Walmod PS
Journal of neuroscience research (2008) 8610: 2201-13. .

Nerve ending "signal" proteins GAP-43, MARCKS, and BASP1.
Mosevitsky MI
International review of cytology (2005) 245: 245-325. .

Natural N-terminal fragments of brain abundant myristoylated protein BASP1.
Zakharov VV, Capony JP, Derancourt J, Kropolova ES, Novitskaya VA, Bogdanova MN, Mosevitsky MI
Biochimica et biophysica acta (2003) 16221: 14-9. .

The BASP1 family of myristoylated proteins abundant in axonal termini. Primary structure analysis and physico-chemical properties.
Mosevitsky MI, Capony JP, Skladchikova GYu , Novitskaya VA, Plekhanov AYU , Zakharov VV
Biochimie (1997) 796: 373-84. .

Access the online factsheet including applicable protocols at <https://sysy.com/product/246003> or scan the QR-code.



FAQ - How should I store my antibody?

Shipping Conditions

- All our antibodies and control proteins / peptides are shipped lyophilized (vacuum freeze-dried) and are stable in this form without loss of quality at ambient temperatures for several weeks.

Storage of Sealed Vials after Delivery

- **Unlabeled** and **biotin-labeled antibodies** and **control proteins** should be stored at 4°C before reconstitution. **They must not be stored in the freezer when still lyophilized!** Temperatures below zero may cause loss of performance.
- **Fluorescence-labeled antibodies** should be reconstituted immediately upon receipt. Long term storage (several months) may lead to aggregation.
- **Control peptides** should be kept at -20°C before reconstitution.

Long Term Storage after Reconstitution (General Considerations)

- The storage freezer must not be of the frost-free variety ("no-frost freezer"). This cycle between freezing and thawing (to reduce frost-build-up), which is exactly what should be avoided. For the same reason, antibody vials should be placed in an area of the freezer that has minimal temperature fluctuations, for instance towards the back rather than on a door shelf.
- Aliquot the antibody and store frozen (-20°C to -80°C). Avoid very small aliquots (below 20 µl) and use the smallest storage vial or tube possible. The smaller the aliquot, the more the stock concentration is affected by evaporation and adsorption of the antibody to the surface of the storage vial or tube. Adsorption of the antibody to the surface leads to a substantial loss of activity.
- The addition of glycerol to a final concentration of 50% lowers the freezing point of your stock and keeps your antibody at -20°C in liquid state. This efficiently avoids freeze and thaw cycles.

Product Specific Hints for Storage

Control proteins / peptides

- Store at -20°C to -80°C.

Monoclonal Antibodies

- **Ascites** and **hybridoma supernatant** should be stored at -20°C up to -80°C. **Prolonged storage at 4°C is not recommended!** Unlike serum, ascites may contain proteases that will degrade the antibodies.
- **Purified IgG** should be stored at -20°C up to -80°C. Adding a carrier protein like BSA will increase long term stability. Many of our antibodies already contain carrier proteins. Please refer to the data-sheet for detailed information.

Polyclonal Antibodies

- **Crude antisera:** With anti-microbials added, they may be stored at 4°C. However, frozen storage (-20°C up to -80°C) is preferable.
- **Affinity purified antibodies:** Less robust than antisera. Storage at -20°C up to -80°C is recommended. Adding a carrier protein like BSA will increase long term stability. Most of our antibodies already contain carrier proteins. Please refer to the data-sheet for detailed information.

Fluorescence-labeled Antibodies

- Store as a liquid with 1 : 1 (v/v) glycerol at -20°C. Protect these antibodies from light exposure.

Avoid repeated freeze-thaw cycles for all antibodies!

FAQ - How should I reconstitute my antibody?

Reconstitution

- All our purified antibodies are lyophilized from PBS. To reconstitute the antibody in PBS, add the amount of deionized water given in the respective datasheet. If higher volumes are preferred, add water as mentioned above and then the desired amount of PBS and a stabilizing carrier protein (e.g. BSA) to a final concentration of 2%. Some of our antibodies already contain albumin. Take this into account when adding more carrier protein. For complete reconstitution, carefully remove the lid. After adding water, briefly vortex the solution. You can spin down the liquid by placing the vial into a 50 ml centrifugation tube filled with paper.
- If desired, add small amounts of azide or thimerosal to prevent microbial growth. This is especially recommended if you want to keep an aliquot a 4°C.
- After reconstitution of fluorescence-labeled antibodies, add 1 : 1 (v/v) glycerol to a final concentration of 50%. This lowers the freezing point of your stock and keeps your antibody in liquid state at -20°C.
- Glycerol may also be added to unlabeled primary antibodies. It is a suitable way to avoid freeze-thaw cycles.
- Please refer to our **tips and hints for subsequent storage** of reconstituted antibodies and control peptides and proteins.